# HYASE™ -10X

20717.04

## EN: Indication for use

Medium for removal of cumulus cells.

## Product Description

HYASE™ is a physiological salt buffer containing hyaluronidase, human serum albumin and penicillin G as an antibacterial agent.

For use after diluting 1:10 with either of the following two alternatives:

1) G-MOPS™ PLUS or G-MOPS™ supplemented with G-MM™ or HSA-solution™ and equilibration at +37°C and ambient atmosphere.

2) G-GAMETE™ and equilibration at +37°C and 6% CO<sub>2</sub> atmosphere.

The concentration of hyaluronidase in HYASE™ is 800 IU/ml before dilution.

## Storage instructions and stability

Store dark at +2 to +8°C.

HYASE™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

## Directions for use

The following is the general procedure for using HYASE™. For more information on recommended use of HYASE™, please see "G Series Manual" or "Vitrolife Fertility Systems Manual". The manuals are available on request from Vitrolife, and may be ordered through the company's internet website, www.vitrolife.com. Direction in other language can also be found on the website.

Dilute HYASE™ either with supplemented G-MOPS™/G-MOPS™ PLUS/G-GAMETE™ to a ratio of 1:10, mix well.

Prepare and label sterile, non toxic multi-well dishes and prepare 3 wells of either supplemented G-MOPS™/G-MOPS™ PLUS/G-GAMETE™ for every well of HYASE™.

Pre-warm the multi well dishes with diluted HYASE™ and the supplemented G-MOPS™/G-MOPS™ PLUS for 15 minutes.

If G-GAMETE™ is used as diluting medium, equilibrate the dishes at +37°C and 6% CO<sub>2</sub> for 3 hours.

Using a sterile, large bore pipette, place 3–5 oocytes into the HYASE™. Gently pipette the HYASE™ and the oocytes to disperse the cumulus cells. Do not expose oocytes to HYASE™ solution for more than 30 seconds.

Use a sterile pipette to transfer oocytes from HYASE™ to the first washing well. Using a fine bore denudation pipette, aspirate each oocyte singly up and down to remove the corona. The diameter of the pipette

should be slightly larger than that of the oocyte. Using a pipette that is too narrow or has a jagged edge may damage the oocyte.

Use the remaining two washing wells to rinse the oocyte. Transfer oocyte to new supplemented G-MOPS™/G-MOPS™ PLUS/G-GAMETE™ or other equivalent droplets for oocyte handling or, if needed, to culture medium for further incubation and maturation.

## Specifications

Sterile filtered	SAL 10 <sup>3</sup>
Bacterial endotoxins (LAL assay) [EU/ml]	< 0.5*
*after diluting 1:10 in G-MOPS™ PLUS	
Cumulus Cell Removal Assay	Positive
Mouse Embryo Assay	Not performed
LOT specific test results are available on the Certificate of Analysis provided with each delivery.	

## Symbols

	Sterilized using aseptic processing techniques
	Temperature limitation
	Do not re-use, discard after procedure.
	Use by - see label.
	Caution: Consult accompanying documents.
	Catalog number
	Batch code

## Precautions

Do not use HYASE™ if it appears cloudy.

HYASE™ contains hyaluronidase, human serum albumin and penicillin G.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBc, HCV, and HTLV I/II and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Vitrolife recommends that media be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

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Vitrolife Sweden AB, Box 9080 (Gustaf Werners gata 2), SE-400 92 Göteborg, Sweden.
Vitrolife Inc., 3601 South Inca Street, Englewood CO 80110, USA.
Vitrolife Pty. Ltd., 55 Flemington Rd, North Melbourne VIC 3051, Australia.
Vitrolife Sweden AB